

Temporary crowns

TRANSLUCENT POLYCARBONATE TEMPORARY CROWNS
DESIGNED TO FIT FOR ALL CLINICAL CASES.



ADVANTAGES

- **Universal shade** for all
- Easy to use and easy to **adjust** thanks to its elasticity
- Perfectly tolerated by periodontal tissue.
- **Occlusal anatomy** very close to natural teeth.
- Resistant to the thermal and chemical variations of the oral cavity.
- Multiple shapes

REFERENCES

Kit of temporary crowns

180 translucent polycarbonate temporary crowns RCP180

Recharges de 5 couronnes :

Upper Right central incisors

L : 7,5 mm	H : 10 mm	N° 14	RCP14
L : 7,7 mm	H : 10,5 mm	N° 13	RCP13
L : 8 mm	H : 11 mm	N° 12	RCP12
L : 8,5 mm	H : 11,5 mm	N° 11	RCP11
L : 9 mm	H : 12,5 mm	N° 10	RCP10
L : 9,5 mm	H : 13 mm	N° 101	RCP101
L : 10 mm	H : 14 mm	N° 100	RCP100

Upper Left central incisors

L : 10 mm	H : 14 mm	N° 103	RCP103
L : 9,5 mm	H : 13 mm	N° 102	RCP102
L : 9 mm	H : 12,5 mm	N° 19	RCP19
L : 8,8 mm	H : 11,5 mm	N° 18	RCP18
L : 8 mm	H : 11 mm	N° 17	RCP17
L : 7,7 mm	H : 10,5 mm	N° 16	RCP16
L : 7,5 mm	H : 10 mm	N° 15	RCP15

Upper Right lateral incisors

L : 5,7 mm	H : 8,5 mm	N° 24	RCP24
L : 6,2 mm	H : 9 mm	N° 23	RCP23
L : 6,5 mm	H : 9,5 mm	N° 22	RCP22
L : 6,8 mm	H : 10 mm	N° 21	RCP21
L : 7mm	H : 10,5mm	N°20	RCP20
L : 7,3 mm	H : 11 mm	N° 2	RCP2

Upper Left lateral incisors

L : 7,3 mm	H : 11 mm	N° 200	RCP200
L : 7 mm	H : 10,5 mm	N° 29	RCP29
L : 6,8 mm	H : 10 mm	N° 28	RCP28
L : 6,5 mm	H : 9,5 mm	N° 27	RCP27
L : 6,2 mm	H : 9 mm	N° 26	RCP26
L : 5,7 mm	H : 8,5 mm	N° 25	RCP25

Upper Right canines

L : 9 mm	H : 13 mm	N° 300	RCP300
L : 8,5 mm	H : 12,5 mm	N° 301	RCP301
L : 8,2 mm	H : 12 mm	N° 30	RCP30
L : 8 mm	H : 11 mm	N° 31	RCP31
L : 7,8 mm	H : 10,5 mm	N° 32	RCP32
L : 7,5 mm	H : 10 mm	N° 33	RCP33
L : 7,2 mm	H : 9,8 mm	N° 34	RCP34

Upper Left canines

L : 7,2 mm	H : 9,8 mm	N° 35	RCP35
L : 7,5 mm	H : 10 mm	N° 36	RCP36
L : 7,8 mm	H : 10,5 mm	N° 37	RCP37
L : 8 mm	H : 11 mm	N° 38	RCP38
L : 8,2 mm	H : 12 mm	N° 39	RCP39
L : 8,5 mm	H : 12,5 mm	N° 302	RCP302
L : 9 mm	H : 13 mm	N° 303	RCP303

Upper premolars

L : 7,4 mm	H : 10 mm	N° 40	RCP40
L : 7 mm	H : 9,5 mm	N° 41	RCP41
L : 6,8 mm	H : 9 mm	N° 42	RCP42
L : 6,5 mm	H : 8,5 mm	N° 43	RCP43
L : 6,2 mm	H : 8 mm	N° 44	RCP44

Lower premolars

L : 7,5 mm	H : 10,2 mm	N° 50	RCP50
L : 7 mm	H : 10 mm	N° 51	RCP51
L : 6,8 mm	H : 9,3 mm	N° 52	RCP52
L : 6,6 mm	H : 9 mm	N° 53	RCP53
L : 6,2 mm	H : 8,5 mm	N° 54	RCP54

Lower incisors (long)

L : 6,2 mm	H : 12 mm	N° 60	RCP60
L : 5,9 mm	H : 11,5 mm	N° 61	RCP61
L : 5,6 mm	H : 11 mm	N° 62	RCP62
L : 5,3 mm	H : 10,5 mm	N° 63	RCP63
L : 5 mm	H : 10 mm	N° 64	RCP64

Lower incisors (short)

L : 5 mm	H : 8 mm	N° 65	RCP65
L : 5,1 mm	H : 8,5 mm	N° 66	RCP66
L : 5,3 mm	H : 9,2 mm	N° 67	RCP67
L : 5,6 mm	H : 10 mm	N° 68	RCP68
L : 6 mm	H : 10,5 mm	N° 69	RCP69